

VELMEX, INC.

7550 Routes 5 & 20 Bloomfield, N. Y. 14469

Telephones: 800 642-6446; 585 657-6151; 585 624-1080

FAX: 585 657-6153 www.velmex.com

Jan. 22. 2008

Converting mounting holes to Metric threaded holes

Product	No.of holes	English Thread Size	Converts to ..
UniSlide A15	4	6/32	M4
UniSlide A25	4	8/32	M5
UniSlide A40	4	10/32	M6
UniSlide A60	4	1/4-20 & 10/32	M8 & M6
UniSlide B90	4	1/4-20 & 5/16-16	M8 & M10
BiSlide	8	1/4-20	M8*
XSlide	4	6/32	M4

*No charge. Spec Metric BiSlide carriage in place of Std.
Typical charge is \$5/hole